



Purchase Quotation

PLEASE REFER TO THIS NUMBER ON
ALL CORRESPONDENCES AND ORDERS

Quote #: Q-404962

Status: Approved

Quote Expiration Date: 10/9/2025

TO:

CUSTOMER NAME	CUSTOMER NUMBER
CONTRA COSTA REGIONAL MEDICAL CENTER	106696
BILL TO ADDRESS	SHIP TO ADDRESS
2500 ALHAMBRA AVE MARTINEZ CA US 94553	2311 LOVERIDGE RD PITTSBURG CA US 94565

TAX INFO:

Hologic is required by law to collect all state and local taxes on all sales. If an exemption certificate is not provided by customer at time of order, final invoices will include these amounts. Many states require both specific operator qualifications and/or licensing and registration of x-ray devices. Hologic is not responsible for fulfilling customer's regulatory obligations.

This Quotation is based on the information known by Hologic regarding your needs and is subject to change or withdrawal by Hologic prior to acceptance. Notwithstanding the foregoing, Hologic may cancel any signed Quotation and/or Customer submitted purchase order if Customer does not take delivery of the Products quoted within one (1) year of execution. This Quotation and the governing terms as noted herein shall supersede all other quotations, agreements, understandings, warranties and representations, whether written or oral, between us, and may be accepted only in accord with their terms. In the event of a conflict between this Quotation and the governing terms, this Quotation shall prevail. To accept, please have an authorized representative sign this Quotation and/or submit a purchase order to your Hologic Representative or to BSH Sales Support as listed below:

BSH Sales Support:

HOLOGIC SALES AND SERVICE, LLC

250 Campus Drive

Marlborough, MA 01752

Fax: (203) 731-8463

BSHSalesSupportUS@hologic.com

ATTN: Angela Womble

Phone: 9253136740

Fax:

Email:

angela.womble@cchealth.org

Quote Date	Requested Date	Hologic Representative	Quote Currency
11/15/2024	1/5/2026	Colleen Morrison colleen.morrison@hologic.com 415-261-7111	U.S. Dollar

Summary of Governing Terms/Contracts	Contract Number	FOB	Payment Terms	Freight Terms
VIZIENT (XR0653) - MAMMO	XR0653	ORIGIN	45 NET	NO CHARGE

3Dimensions System - Pittsburg

Qty	Product Name	Description	List Price	Unit Price	Extended Price
1	3DM-SYS-STD	3DIMENSIONS MAMMOGRAPHY SYSTEM 3D US	\$843,000.00	\$403,000.00	\$403,000.00
1	PRD-04420	HIGH RESOLUTION READY DETECTOR	Included	Included	Included
1	3DM-KIT-ERGO-XRAY	3DM UAWS ERGONOMIC UPGRADE	\$11,800.00	\$5,150.00	\$5,150.00
1	DIM-LIC-I2D3DQ	INTELLIGENT 2D AND 3DQUORUM LICENSES FOR NEW 3DIMENSIONS SYSTEM	\$92,900.00	\$51,000.00	\$51,000.00
1	3DQ-TRAIN-INIT-01	3DQUORUM, INITIAL TRAINING, VIDEO-ON-DEMAND, UNLIMITED TECHNOLOGISTS AND PHYSICIANS	Included	Included	Included
1	DIM-LIC-GDT-IC-V2	GENIUS AI DETECTION 2.0 AND IMAGECHECKER CAD SOFTWARE LICENSES	\$77,900.00	\$42,000.00	\$42,000.00

Qty	Product Name	Description	List Price	Unit Price	Extended Price
1	GDT-TRAIN-INIT-01	GENIUS AI DETECTION, INITIAL TRAINING, ON-DEMAND, UNLIMITED TECHNOLOGISTS AND PHYSICIANS	Included	Included	Included
3	ASY-04662	RACK, PADDLE STORAGE	\$471.00	\$333.00	\$999.00
1	PHANTOM-ACR-156	ACR 156 PHANTOM	\$1,220.00	\$922.50	\$922.50
1	SVC-SDM-OPT-BTO	CONFIGURE SYSTEM OUTPUT TO BTO FORMAT	Included	Included	Included
1	DIM-TRAIN-APPS-INIT	TECHNOLOGISTS, DIMENSIONS, INITIAL TRAINING, 2 DAYS, 1 SITE, MAX 5 TECHNOLOGISTS	Included	Included	Included
1	DIM-TRAIN-PHY-INIT	MEDICAL PHYSICIST, DIMENSIONS, INITIAL TRAINING, 8 HRS (5 HRS LIVE - 3 HRS ONLINE TRAINING), 1 SITE, MAX 2 PHYSICISTS	Included	Included	Included
1	DIM-TRAIN-RAD-INIT	RADIOLOGISTS, TOMOSYNTHESIS, INITIAL TRAINING, 8 HOURS VIRTUAL TRAINING, 14 RADIOLOGISTS	Included	Included	Included
1	TRADE-REC-GE_ESS2D	TRADE-IN OF GE SENOGAPHE ESSENTIAL 2D SYSTEM (INSTALLED AND FUNCTIONING)	Included	\$-25,000.00	-\$25,000.00
1	HLX-COMPETITIVE-DISCOUNT	HOLOGIC COMPETITIVE DISCOUNT. Please see the Product Long Description for terms and conditions.	Included	\$-20,000.00	-\$20,000.00
1	SCRAP-TIER1-FEE-BSH	SCRAPPING FEE IS APPLICABLE TO TIER1 BSH CAPITAL ITEMS	\$2,500.00	\$2,500.00	\$2,500.00
3Dimensions System - Pittsburg TOTAL:					\$460,571.50

Portfolio Volume Discount

Qty	Product Name	Description	List Price	Unit Price	Extended Price
1	BSH-VOLUME-DISCOUNT	BSH PORTFOLIO VOLUME DISCOUNT	Included	\$-10,000.00	-\$10,000.00
Portfolio Volume Discount TOTAL:					\$-10,000.00

List Price Total:	USD 1,030,733.00
Discount:	USD 580,161.50
Total Quote Price:	USD 450,571.50
Tax:	USD 35,156.62
Final Quote Price:	USD 485,728.12

Customer agrees to keep the discount price provided to them in this Quotation or agreement confidential and not disclose it to anyone other than as required by law or court order.

The BSH Portfolio Volume Discount for multi-product purchases is offered to Customer for purchase of the Products herein. Customer is prohibited from returning any Portfolio Volume Discount products purchased hereunder except as outlined in the applicable terms and conditions.

Any disputes arising out of the purchase of the Products quoted herein will be governed by the laws of the state of California. In the event of a conflict with this note and any other applicable terms, this note will control.

Title Transfer

The undersigned CONTRA COSTA COUNTY ("Transferor"), for valuable consideration hereby transfer title for property defined as:

3Dimensions System - Pittsburg

Product Name	Trade-In Manufacturer	Trade-In Serial #
TRADE-REC-GE_ESS2D	GE	606736BU7

to Hologic Sales and Service, LLC of 250 Campus Drive, Marlborough, MA 01752, effective as of date of the actual pick-up of the property by Hologic Sales and Service, LLC.

If applicable, Transferor acknowledges as such consideration a trade-in credit value of

US \$ 25,000.00

Transferor warrants in regard to the property being transferred that:

1. Transferor is the legal owner,
2. Transferor will defend the title against any and all claims and demands of all persons,
3. Transferor will take all further steps necessary to effectively transfer ownership to Hologic Sales and Service, LLC, at no cost or expense to Hologic Sales and Service, LLC, and
4. It is free from all liens and encumbrances
5. All patient data stored within the hard drive(s) of transferred property has been removed from all data storage devices in accordance with the Health Insurance Portability and Accountability Act of 1996 ("HIPAA") and the Health Information Technology for Economic and Clinical Health of 2009 ("HITECH") regulations implemented by the U.S Dept. of Health and Human Services.
6. The transferred property has been thoroughly cleaned and, if property has been in contact with potentially infectious materials (blood, blood products or other potentially infectious materials), it must be decontaminated in accordance with OSHA 29 CFR 1910.1030 "Bloodborne Pathogens." If it is not feasible to fully decontaminate the equipment, the equipment must be clearly labeled with a warning sign that indicates the possible presence of blood/body fluids. Decontamination must be done with an appropriate EPA approved disinfectant - <https://www.epa.gov/pesticide-registration/selected-epa-registered-disinfectants>).

If Transferor and Hologic are parties to a Service Agreement providing maintenance and repair services for the property being transferred herein, then this Title Transfer shall serve as an amendment to said Service Agreement terminating coverage with respect to the transferred property. The termination of coverage shall be effective as of the date of the actual pick-up of the property by Hologic Sales and Service, LLC.

To HAVE AND TO HOLD the same, with the rights, privileges, and appurtenances thereof, unto Hologic Sales and Service, LLC, its successors and assigns, forever, to its own use and behalf.

Upon receipt of a purchase order and/or signed Quote, your Hologic team will work collaboratively on an installation timeline.

Sales Orders that are requested to be cancelled within forty-five (45) days of the confirmed installation date must be approved by Hologic and may be subject to a cancellation fee of ten percent (10%) of the total Quote price for the items contained herein.

Once the installation confirmation is provided by Hologic, all requests to reschedule an installation within seven (7) business days of the confirmed installation date may be subject to a rescheduling fee of \$2,500.00 USD.

Buyer Acceptance

CONTRA COSTA COUNTY

By: _____ (signature)

Name: _____ Title: _____ (print/type)

Date: _____

_____ I will issue a PO for this purchase. Please do not ship the Products until a PO is issued.

_____ I will not issue a PO for this purchase. A PO is not required for shipment and invoicing. Please ship the Products based on this signed Quotation.

Additional Buyer Acceptance (if applicable)

By: _____ (signature)

Name and Title: _____ (print/type)

Date: _____

Please provide the Shipping and Billing address here if different from the quote address above
(If this section is left blank, the product will ship and bill to the addresses printed at the top)

Shipping Address

Billing Address:

Hologic Approval:



Date:

Product Name	Long Description
3DM-SYS-STD	3Dimensions™ mammography system for Genius®3D Mammography screening and diagnostic mammography with 2D capabilities. Includes Hologic Clarity HD™ technology and SmartCurve™ breast stabilization system. Upgradeable to interventional imaging. INCLUDES: X-ray Gantry: • Generator: Fully integrated constant potential, high frequency, inverter type. • Detector: Hologic Clarity HD High-resolution enabled detector. • X-ray Tube: Tungsten, bi-angular, high speed, high heat capacity. • X-ray Filters: Rhodium, silver, aluminum. • Anti-scatter Grid: Auto-retracting linear grid. • Dual-function gantry footswitches (2).. Adjustable-height Acquisition Workstation: • CPU: High-performance computer, multi-core Intel-based CPU, minimum 32 GB RAM, minimum 4 TB disk, Windows 10/64, high-performance NVIDIA GPU. • Includes DVD +/- R/W. • User Interface: 1.2 MP color LCD touchscreen. • Image Review: 21.3" 3 MP, medical-grade color, DICOM monitor, double-jointed articulating arm. • Full X-ray shield, X-ray exposure footswitch, pull-out keyboard drawer, keyboard and mouse, Uninterruptible Power Supply (UPS). 3Dimensions System Software: • User access control, patient and study selection, imaging procedure selection and definition, X-ray parameter control, image review and acceptance/rejection, quality control. • Licenses: Hologic Clarity HD Imaging, Tomosynthesis Imaging, SmartCurve System, Diagnostic Imaging, Dynamic Tube Head Motion, Advanced Connectivity, Bi-directional Communication (Notices). Connectivity: • DICOM: Modality work list, storage, storage commitment, query/retrieve, print. • IHE Profiles: Scheduled workflow, patient information reconciliation, mammography image. • Advanced Connectivity: Ability for the system to participate in DICOM Modality Performed Procedure Step (MPPS) transactions and to output DICOM Radiation Dose Structured Report (RDSR) objects to third-party dose aggregation and reporting systems. The MPPS and RDSR capabilities can be enabled and disabled independently. Accessories: • SmartCurve System (3): Small 18 x 24 cm, large 24 x 29 cm, Mini 18 x 4. • Screening Paddles (3): 18 x 24 cm, 24 x 29 cm, small breast. • Diagnostic Paddles (5): Frameless spot contact, 7.5 cm spot contact, 10 cm spot contact, 10 cm magnification, 7.5 cm spot magnification. • Magnification stand with platform, flat field phantom and case, ACR tomosynthesis geometry calibration phantom, Dimensions interconnect cable kit, 2D fixed face shield, retractable 3DMammography face shield. User, service, maintenance, and QC manuals. MammoPad® Breast Cushion trial box Installation: Installation by Hologic certified technicians. Warranty: Standard one-year parts and labor warranty. Warranty service coverage is Monday-Friday from 8:00 am to 5:00 pm local time. CONDITIONS: • Required but independently sold software features are not included gratis under Hologic Warranty, or Hologic Service Contracts that include software upgrades, and must be ordered separately. • Required system software updates that are included under Hologic Warranty, or Hologic Service Contracts will be provided at no charge.

Product Name	Long Description
PRD-04420	Hologic 3D Mammography™ high-resolution ready digital image receptor for 3Dimensions™ mammography systems or Selenia® Dimensions® systems. INCLUDES: • Digital Image Receptor • Amorphous selenium, TFT • Structure: Single 24 x 29 cm plate • Image Matrix Sizes: 2560 x 3328 (18 x 24 cm); 3328 x 4096 (24 x 29 cm) •Pixel Size: 0.070 mm • Limiting Spatial Resolution: 7.1 lp/mm. (Hologic Clarity HD license purchase required to enable Clarity HD high-resolution imaging on 3D Performance and Dimensions 6000 and 9000 systems. This license is included with a 3Dimensions system.)
3DM-KIT-ERGO-XRAY	Provides an integrated ergonomic package for the 3Dimensions™ system. Includes: • Fingerprint reader for quick biometric login of the technologist, manager or other site personnel • Powered memory height adjust • Bar code scanner Note: • Kit only available for 3Dimensions systems

Product Name	Long Description
DIM-LIC-I2D3DQ	3DQuorum® and Intelligent 2D® technology for new 3Dimensions™ systems. 3DQuorum technology utilizes Genius AI®-powered analytics to uniquely reconstruct 1mm high-resolution 3D™ data to produce 6mm high-resolution SmartSlices. SmartSlices are designed to expedite reading time by reducing the number of images to review. Intelligent 2D technology generates 70µm resolution synthesized 2D images that are smart mapped to the 3D™ data. INCLUDES: • 3DQuorum and Intelligent 2D software licenses. REQUIREMENTS: • Hologic Clarity HD® high resolution 3D™ imaging technology. • Tomosynthesis license. • High-performance AWS computer (provided if needed). • Dimensions software minimum version 1.10/2.1 with Windows 10 OS. • BTO image data format is required for 3DQuorum technology. 3DQuorum software is unable to process SCO data. NOTES: • Hologic SecurView® workstation software 10.4 minimum. • For customers reading on PACS workstations minimum requirements will vary depending on the vendor. • Licenses are non-transferable to other mammography systems. • The licenses do not extend the Dimensions system warranty. New hardware, if needed, has a replacement part warranty as per Hologic Sales Terms and Conditions.
3DQ-TRAIN-INIT-01	3DQuorum™ imaging technology training provided via Video-On-Demand from www.hologic.com/training and/or for physicians from HologicAce.com. The Hologic Ace portal also has a DICOM viewer and images for physicians to review. Training may be accessed at any time for unlimited staff for an unlimited period. Initial training is included in the purchase price of your license.

Product Name	Long Description
DIM-LIC-GDT-IC-V2	The Genius AI™ Detection 2.0 software applies deep learning AI to tomosynthesis images from ClarityHD™ and standard resolution images. Results can be displayed on tomosynthesis slices, 3DQuorum™ SmartSlices, synthesized 2D images (C-View™ and Intelligent 2D™), and FFDM images (Only when associated with Combo or ComboHD imaging). The ImageChecker® CAD software applies markings to standalone FFDM cases only. For new gantry purchases. INCLUDES: • Genius AI Detection 2.0 software license • Indicates suspicious areas with lesion score • Prioritizes most suspicious cases • Complexity score mode batches reading cases according to complexity • Read time indicator can distribute cases among radiologists • Real time case results available at the acquisition workstation • Correlate lesions between CC and MLO views (Compatibility varies by PACS vendor) • ImageChecker CAD software license • Citra advanced CAD display license provides additional information about why ImageChecker CAD marked specific regions: • RightOn™ CAD marks indicate the types of features that were detected • PeerView™ CAD marks provide anatomic outline of tissues • EmphaSize™ markers are scaled according to the prominence of features • Installation on one (1) new Selenia® Dimensions® or 3Dimensions™ AWS • On-demand webinar training available for Genius AI Detection 2.0 • In-person training, upon request, available in conjunction with gantry installation on SecurView REQUIREMENTS: • License installed with new gantry • AWS software minimum 1.12/2.3 (on minimum computer CMP-01816 with Windows 10 software) • Hologic workstation software minimum: SecurView® workstation 11.0 • Hologic workstation software recommended: SecurView® workstation 12.0 NOTE: • Hologic works closely with the major PACS companies to provide sample test images with corresponding results and guidelines on implementing the most efficient workflow for Genius AI Detection software. We recommend that you speak directly with your PACS vendor to determine how Genius AI Detection results will be displayed in your PACS environment. If your PACS vendor needs additional information, please provide the PACS vendor contact information to your Hologic Account Executive. • Advanced Citra CAD features require workstations that conform to proper display of those features. They can be disabled for use with other non-conformant workstations. The customer needs to check with their workstation vendor.
GDT-TRAIN-INIT-01	Initial training is provided On-Demand from www.hologic.com/training and/or for physicians from HologicAce.com. Training may be accessed at any time for unlimited staff for an unlimited period. Initial training is included in the purchase price of your license.

Product Name	Long Description
ASY-04662	Wall-mounted, felt-lined rack for storing system compression paddles. Each rack, sold separately, provides enough storage room for 3-4 compression paddles. Includes: • Paddle storage rack with felt lining • Wall-mounting bracket (installation not included) • Graphic paddle labels Dimensions: • W x H x D: 36 1/4" x 7" x 4" (from the wall) • Minimum Vertical Pitch: 12" to 14" when multiple racks are installed one above another Recommended: • Selenia® Dimension® Avia systems: min. 1 rack • Selenia Dimensions 2D systems: min. 2 racks • Selenia Dimensions 3D™ systems: min. 3 racks • 3Dimensions™ systems: min. 3 racks Requirements: • Must be securely attached to the wall • Must be installed by a professional installer
PHANTOM-ACR-156	The Mammographic Accreditation Phantom manufactured by Gammex is designed to test the performance of a mammography system's image quality and sensitivity using target objects in the phantom to simulate calcifications, fibrous calcifications in ducts, and tumor masses. The phantom simulates the X-ray attenuation of a 4.2 cm compressed human breast composed of 50% adipose tissue and 50% glandular tissue. Target objects within the phantom range in size, shape, and density, similar to those found clinically. • Breast phantom is compatible with digital and analog equipment. • Approved by ACR for Mammography. Image quality and system sensitivity follow ACR and MQSA guidelines. • Dimensions: Height 1.75 in. (4.5 cm) x width 4 in. (10.2 cm) x depth 4.25 in. (10.8 cm)
SVC-SDM-OPT-BTO	This configuration enables output of tomosynthesis slices in DICOM Breast Tomosynthesis Image Object form. Use of this configuration will first require an integrated planning team, including your IT department, Hologic and other vendors, to work together to ensure that your enterprise is ready for use of the tomosynthesis data in DICOM Breast Tomosynthesis Image Object form. The completion of critical feasibility questions included in Hologic's Enterprise Survey will guide the team through understanding any infrastructure requirements and changes necessary. Note: While a preliminary check by your Hologic representatives may have allowed the ability to quote this output configuration, completion of the Enterprise Survey is required before Dimensions Tomosynthesis system or option is enabled. Software and hardware upgrades may be required. Hologic makes no guarantees of software and hardware performance for products not associated to Hologic. By signing this quote, the customer agrees that the completion of purchase of the accompanying Hologic products shall not be contingent on the implementation of this no-charge configuration. Requires: • Completion of Hologic Enterprise Survey by site personnel in conjunction with Hologic representatives • PACS system including Deep Archive capable of storing / retrieving DICOM Breast Tomosynthesis Image Objects and with suitable storage capacity • Softcopy review workstation capable of displaying DICOM Breast Tomosynthesis Image Objects

Product Name	Long Description
DIM-TRAIN-APPS-INIT	Technologist training for a new 3Dimensions™ or Selenia® Dimensions® system. Initial training is included in the purchase price of your system. Training duration to be determined by customer needs up to 2 days. • One session of technologist onsite applications training or other clinical support for maximum of 5 technologists based on training effectiveness and space limitations. • Additional sessions may be needed. For additional groups of up to 5 technologists add purchasable DIM-TRAIN-APPS-ADDL. • Online CEU courses required prior to onsite training. • Video training available during and post training. • Onsite portion of training or other clinical support must be completed within 24 months of equipment installation. Required FDA training: • FFDM accreditation is required: apply to the ACR or your State for FFDM certification. • Once FFDM accredited, contact the MQSA FFDM Certification Extension Program for Tomosynthesis. • Sites must obtain ACR or State FFDM accreditation before using the tomosynthesis modality. Please note: Cancellation must be made 6 business days prior to the confirmed and scheduled training start date. Fee for cancellation within 6 business days of confirmed training is \$2,000.
DIM-TRAIN-PHY-INIT	Medical physicists training for a 3Dimensions™ or Selenia® Dimensions® system. Initial training is included in the purchase price of your system and is valued at \$1,500 (unused training cannot be deducted from your purchase price). Five (5) hours of live tomosynthesis training for up to 2 medical physicists with a Field Service Engineer during the installation of the system and access to Hologic's 3-hour online tomosynthesis training course for medical physicists to fulfill the 8-hour FDA requirement. See course description for more details and requirements. Visit www.hologic.com/training for a complete list of Hologic educational opportunities. Please note: Cancellation must be made 6 business days prior to the confirmed and scheduled training start date. Fee for cancellation within 6 business days of confirmed training is \$2,000.
DIM-TRAIN-RAD-INIT	Radiologist training for new 3Dimensions™ or Selenia® Dimensions® system(s). Initial training is included in the purchase price of your system and is valued at \$5,250 (unused training cannot be deducted from your purchase price). Access to Hologic's virtual tomosynthesis training course for up to 14 radiologists. This program fulfills the 8-hour FDA training requirement. See course description for more details and requirements. Visit www.hologic.com/training for a complete list of Hologic educational opportunities.
HLX-COMPETITIVE-DISCOUNT	Customers who trade in their existing competitive mammography gantry for a new Hologic 3D™ capable system, are eligible for an additional "Competitive Discount" in conjunction with their trade-in credit. Please see notes section for additional eligibility requirements.
SCRAP-TIER1-FEE-BSH	SCRAPPING FEE IS APPLICABLE TO TIER1 BSH CAPITAL ITEMS FOR ALL RELATED COSTS INCLUDING BUT NOT LIMITED TO SCRAP FEE, DISPOSAL FEE, HANDLING FEE AND OTHERS.



California Radiation Control Regulations Requirements Form

NOTICE OF REQUIREMENT TO REGISTER FACILITY & EQUIPMENT

Title 17, Section 30118(4)(b) of the California Code of Regulations requires “the vendor shall inform the receiver of each machine of the registration requirements of Section 30108 of these regulations”. Therefore, Hologic, Inc. is notifying you of the following:

Facility & Equipment Registration:

California state law requires facilities & equipment to be registered within thirty (30) days of acquisition of x-ray equipment with the state Department of Health Services. This includes modifications to the facility and/or addition of new equipment. (See California Code of Regulations, Title 17, Sections 30108 through 30146).

Certain counties require notification and/or inspections prior to installation. You may also have to register with your local Dept. of Health Services.

Bone Density Equipment (Operators)

A radiologist with a Radiology Supervisor license or a physician with a Radiography Supervisor license, a CRT, or a RT can operate the DXA machine. Non-licensed and non-RT operators are required to obtain a Limited Permit to operate an X-ray Bone Densitometer. Please call the California Department of Health Services, Phone: (916)327-5106, to obtain an x-ray bone densitometry limited permit application form. To obtain a Limited Permit operators are required to attend a training program approved by the State of CA, followed by an exam.

It is your responsibility to contact the appropriate state and local agencies to obtain the required forms and instructions. Please see attached page for contact information.

Acknowledgement of notice

Please complete the following information and return to your Hologic, Inc. contact.

Facility Name:			
Street Address:			
City:		Zip Code:	
Name:			
Title:		Phone Number:	
Signature:		Date:	

California Radiation Control Regulations Requirements Form

California

<https://www.cdph.ca.gov/Programs/CEH/DRSEM/Pages/RHB.aspx>

Office	Address	Telephone/Fax Numbers
Radiologic Health Branch Sacramento	Mailing Address: Department of Health Services Radiologic Health Branch PO Box 997414, MS 7610 Sacramento, CA 95899-7414 Physical Address: Department of Health Services Radiologic Health Branch 1500 Capitol Avenue, 5 th Floor, MS 7610 Sacramento, CA 95814-5006	Phone: (916) 327-5106 Fax: (916) 440-7999
Radiologic Health Branch Richmond	850 Marina Bay Parkway, Bldg. P, 1 st Floor Richmond, CA 94804	Phone: (510) 620-3416 Fax: (510) 620-3874
Radiologic Health Branch Brea	500 S. Kraemer Blvd Brea, CA 92821	
	Radioactive Materials (RAM), Suite 235	Phone: (714) 524-1409 Fax: (714) 524-1908
	X-ray (Region 2*), Suite 235	Phone: (714) 524-4450 Fax: (714) 524-1908
	X-ray (Region 5*), Suite 225	Phone: (714) 524-5681 Fax: (714) 524-5682
Radiologic Health Branch Log Angeles County	Radiation Management 3530 Wilshire Boulevard, 9 th Floor Los Angeles, CA 90010	Phone: (213) 351-7897 Fax: (213) 351-2718
Radiologic Health Branch San Diego County	Mailing Address: County of San Diego Department of Environmental Health, Radiological Health Program 5500 Overland Ave, Suite 110 MS O-560 San Diego, CA 92123	Phone: (858) 694-3621 Fax: (858) 694-3629

* **Region 2** serves the following counties:

Inyo, Kern, Kings, Mono, Monterey, San Benito, San Luis Obispo, Santa Barbara, Santa Cruz, Tulare, Ventura.

* **Region 5** serves the following counties:

Imperial, Orange, Riverside, San Bernardino.